

CE 2460



HANGZHOU TECHDERM BIOLOGICAL PRODUCTS CO., LTD

Floor 2, East Building 6, No.88 Jiangling Road,
Binjiang District, Hangzhou, P.R.China

Tel: 86 571 87165526

Fax: 86 571 87165526

Email: info@techderm.com

Website: www.techderm.com



Landlink GmbH

ADD: Dorfstrasse 2/4, 79312, Emmendingen, Germany

TEL: 0049 7641 9626855

FAX: 0049 7641 9626855



PROFESSIONAL MANUFACTURER OF
HYALURONIC ACID GEL

Sofiderm NATURE BEAUTY JUST STARTS HERE

Sofiderm NATURE BEAUTY JUST STARTS HERE

Sofiderm



1. What is Sofiderm

Sofiderm is a kind of hyaluronic acid gel for plastic surgery, which can be used for smoothing away all kinds of wrinkles and body shaping. It is non-animal ingredient and has 6-24 months' longevity.




2. How Does Sofiderm Work?

Once injected into the skin, Sofiderm gel will block skin space and restore lines and wrinkles quickly. Skin rejuvenation will start easily and end quickly. This easy operation will just take a few minutes.

MCLPE (Multiphase Cross Linking & Purification Enhancement) Technology

MCLPE technology is an innovative and advanced technology which combines high cross-linking efficiency and high degree of purity of natural HA, resulting in a 3D gel with firmer texture and long-lasting effect, and creates smooth and natural look with a proven safety record.



Unique particle size to meet different requirements



Sofiderm Finelines

Sofiderm Derm

Sofiderm Deep

Sofiderm Derm Plus

Sofiderm Derm Sub-Skin



3. Advantages

ENHANCED SAFETY

The Origin of HA is microbial fermentation and animal free, which ensures its safety. The hyaluronic acid raw material is approved by American FDA and European EDQM, which is reliable and trustworthy.

LONG LASTING

Sofiderm multiterm cross linking technology ensures strong stability of hyaluronic acid gel which brings a long lasting effect.

LOW RATE OF ALLERGIC REACTION

Strict quality control system brings very low rate of allergic reactions.

HIGH SATISFACTION

Patients and doctors are satisfied with the natural results right after the treatment.

Sofiderm Products Range



SOFIDERM FINELINES	
Ingredient :	Cross-linked hyaluronic acid gel
Concentration :	20±3mg/ml
Particle Size :	0.1-0.15mm
Number of particles /ml:	200,000
Volume of syringe:	1ml/2ml
Needle size:	27GX2
Duration:	6-18 months
Injection Depth:	Upper part of dermis
Indications:	Crow's feet④ Corrects scars③⑤ Frown②



SOFIDERM DERM	
Ingredient :	Cross-linked hyaluronic acid gel
Concentration :	20±3mg/ml
Particle Size :	0.15-0.28mm
Number of particles /ml:	100,000
Volume of syringe:	1ml/2ml
Needle size:	27GX2
Duration:	6-18 months
Injection Depth:	Middle part of dermis
Indications:	Lip⑥⑦⑧ Cheek⑪ Marionette lines⑨ Forehead lines①



SOFIDERM DEEP	
Ingredient :	Cross-linked hyaluronic acid gel
Concentration :	20±3mg/ml
Particle Size :	0.28-0.5mm
Number of particles /ml:	10,000
Volume of syringe:	1ml/2ml
Needle size:	27GX1+26GX1
Duration:	6-24 months
Injection Depth:	Deep layer of dermis
Indications:	Chin ⑩ Nasolabial folds ⑫ Temples ⑭ Nose ⑬



Sofiderm Products Range



SOFIDERM DERM PLUS

Ingredient :	Cross-linked hyaluronic acid gel
Concentration :	20±3mg/ml
Particle Size :	0.5-1.25mm
Number of particles /ml:	5,000
Volume of syringe:	10ml/20ml
Needle size:	No needle
Duration:	6-24 months
Injection Depth:	Deep layer of dermis or subcutaneous tissue shallow
Indications:	Breast Enhancement ⑨



SOFIDERM DERM SUB-SKIN

Ingredient :	Cross-linked hyaluronic acid gel
Concentration :	20±3mg/ml
Particle Size :	1.25-2.0mm
Number of particles /ml:	2,000
Volume of syringe:	10ml/20ml
Needle size:	No needle
Duration:	6-24 months
Injection Depth:	Deep layer of dermis or subcutaneous tissue shallow
Indications:	Buttock Enhancement ⑩






SOFIDERM MESO

Ingredient:	Non-crosslinked hyaluronic acid gel
Concentration :	15mg/ml
Volume of syringe:	1ml/2ml/3ml
Needle size :	30GX2
Duration:	1 month
Injection Depth:	Upper layer of dermis
Function:	Skin rejuvenation, skin moisturizing, skin enlightening

LONG LASTING EFFECT


According to the clinical evaluation data, in the 24-week long clinical research of single blind and paired design, the percent of obvious improvement on injection sites (which are nasolabial lines, fishtail wrinkles, forehead wrinkles, depression scar, glabella wrinkles, etc) is from 65.1% to 82.09%.



✿ Results may last up to 18 months after injection with correction still noticeable (60%) after 24 months.

HIGH SATISFACTION

Sofiderm shows high satisfaction from the patients and doctors with the natural volume right after the treatment.



HANGZHOU TECHDERM BIOLOGICAL PRODUCTS CO., LTD

ABOUT OUR FACTORY:

Techderm is a professional HA dermal filler manufacturer in China, which is established in the year of 2012 and located in Hangzhou City, China. We offer top quality cross-linked HA fillers for plastic surgery all over the world. The perfect product development, production, sales and after-sale service system make us outstanding.

ABOUT OUR EXPERIENCE:

Our R&D team has been studying in the field of hyaluronic acid gel for more than 7 years and the team leader, Mr. Chai, has more than 10 years' experience on it. Experienced researchers and workers are the firm guarantee of Hangzhou Techderm.

ABOUT OUR EQUIPMENT:

- 01 GMP Standard Work Shop
- 02 Sterilized Laboratory
- 03 Complete Precision Instruments

